

### INFORMATION

Customer Name: ..... Date: .....  
 Address: ..... PO#: .....  
 City: ..... Phone: .....  
 Account: .....  
 Authorizer Name: .....

### STANDARD NO CHARGE SPECIFICATIONS

SEAT WIDTH: 14" – 18"	SEAT DEPTH [ADJUSTABLE]: 16" – 18"	SEAT HEIGHT [ADJUSTABLE]: 10.25" – 18.25"
BACK ANGLE [ADJUSTABLE]: 0°, 5°, 10°, 15° [DEGREES]	WEIGHT OF CHAIR: 60 LBS [APPROXIMATELY]	SINGLE RELEASE, DUAL CYLINDERS
ADJUSTABLE STROLLER HANDLE: 0° – 180° [DEGREES]	DYNAMIC TILT ADJUSTMENTS 0 TO 25° [DEGREES]	STANDARD LOCKOUT COLLARS
WEIGHT CAPACITY: 250 LBS	TWO-WAY LOCKING GAS STRUT	STANDARD T STYLE ARMS
REAR WHEEL WITH HANDRIMS AND FRONT CASTERS	REAR HEIGHT ADJUSTABLE ANTI TIPPERS	

#### FRAME TYPE O1

- EXTB Dynamic Frame - Black \$2,699
- EXTINTS Independant Power Tilt System (A1) \$2,000

#### CLIENT WEIGHT O2

- EXTWA Up to 250Lbs N/C
- EXTWAMF Exceeds 250Lbs \$730
- EXTWAMG Exceeds 350Lbs \$900
- EXTWAMF Exceeds 450Lbs \$1,110

#### SEAT WIDTH O3

- EXTSW14 Seat Width 14" N/C
- EXTSW15 Seat Width 15" N/C
- EXTSW16 Seat Width 16" N/C
- EXTSW17 Seat Width 17" N/C
- EXTSW18 Seat Width 18" N/C
- EXTSW19 Seat Width 19" (A5) - W111 \$200
- EXTSW20 Seat Width 20" (A5) - W111 \$200
- EXTSW21-24 Specify: ..... (A5) - W111 \$450
- EXTSW25-32 Specify: ..... \$450

#### FOOT REST & FOOT PLATES O4

- EXTOFST Offset Foot Rest N/C
- EXTF6A 60 Degree Hanger  Pivot  Pin N/C
- EXTF7A 70 Degree Hanger  Pivot  Pin N/C
- EXTF8A 80 Degree Hanger  Pivot  Pin N/C
- EXTEA Elevating Legrest  Pivot  Pin (A2) - W106 \$404
- EXTAL Articulating Legrest  Pivot  Pin (A2) - W106 \$450
- EXTNR Non-removable Pin Style Legrest N/C
- EXTEB Elevating Legrest System \$450
- EXTPPL Plastic Footplates N/C
- EXTPXL Extra Large Footplates \$150
- EXTPAF Angle Adjustable Footplates - W107 \$150
- EXT1FB One Piece Foot Board Padded \$80
- EXTFRE Foot Rest Extensions  Short  Medium  Long N/C

#### SEAT DEPTH O5

- EXTSD16 Seat Depth 16" N/C
- EXTSD17 Seat Depth 17" N/C
- EXTSD18 Seat Depth 18" N/C
- EXTSD19 Seat Depth 19" \$200
- EXTSD20 Seat Depth 20" \$200
- EXTSD21-24 Specify: ..... \$450

#### BACK TYPE O6

- EXTAB1 Angle Adjustable Back - Preset N/C
  - 0°  5°  10°  15°
- EXTENDO Endomorph Canes W ( ) N/C
- EXTRBK Recline Kit \$560

#### SEAT TO FLOOR SPECIFICATIONS (± 0.5")

##### (LOW SEAT SETTING)

STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/> 10.50"	16"	3"	2 HOLE
<input type="checkbox"/> 13.00"	20"	5"	2 HOLE
<input type="checkbox"/> 13.75"	20"	5"	2 HOLE
<input type="checkbox"/> 14.50"	20"	6"	3 HOLE
<input type="checkbox"/> 14.50"	22"	5"	3 HOLE
<input type="checkbox"/> 14.50"	22"	6"	3 HOLE
<input type="checkbox"/> 15.25"	20"	5"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 15.25"	20"	6"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 15.25"	22"	5"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 15.25"	22"	6"	<input type="checkbox"/> 3 HOLE <input type="checkbox"/> 5 HOLE
<input type="checkbox"/> 16.00"	22"	6"	5 HOLE

##### (HIGH SEAT SETTING)

STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/> 16.75"	22"	6"	5 HOLE
<input type="checkbox"/> 17.50"	22"	6"	5 HOLE
<input type="checkbox"/> 17.50"	24"	6"	5 HOLE
<input type="checkbox"/> 18.25"	24"	8"	5 HOLE

### ARM REST Select one below 07

- EXTTA T-Arm Ht. Adj./ Removable (A3) N/C
- EXTFA Flip Back Ht. Adj./ Removable N/C
- EXTTAS T-Arm Ht. Adj./ Removable (A3) N/C
- EXTFAS Flip Back Ht. Adj./ Removable N/C
- EXTGWAFP Gel Waterfall Arm Pads (A1) \$150

### REAR WHEELS Select one below 08

- EXTRWTO 20" Solid Treaded Tire N/C
- EXTRWT2 22" Solid Treaded Tire N/C
- EXTRWT4 24" Solid Treaded Tire N/C
- EXTRWT12I 12" Solid Tire Inside Mount 120
- EXTRWT16 16" Solid Tire N/C
- EXTAQR Axles - Quick Release \$120
- EXTHRT Narrow Handrim Config. (1/2" less) N/C
- EXTPCHR Plastic Coated Hand Rims N/C

### CASTER Select one below 09

- EXTFC1D 3" X 1" Caster N/C
- EXTFC1A 5" X 1.25" Caster N/C
- EXTFC1B 6" X 1.25" Caster N/C
- EXTFC1C 8" X 1.25" Caster N/C
- EXTFC2A 5" X 1.75" Caster (Soft Roll) \$80
- EXTFC2B 6" X 1.75" Caster \$80
- EXTFC2C 8" X 1.75" Caster \$80
- EXTFC4B 6" X 2 Caster (Steel Fork) \$80
- EXTFC4C 8" X 2 Caster (Steel Fork) \$80

### WHEEL LOCK Select one below 10

- EXTALPH Aluminium Push To Lock N/C
- EXTALPL Aluminium Pull To Lock N/C
- EXTEXTH Steel Push To Lock N/C
- EXTEXTL Steel Pull To Lock N/C
- EXTBRXC 4" Brake Extension \$40
- EXTBRXA 6" Brake Extension \$40
- EXTBRXB 9" Brake Extension \$130

### SEATING BENIFITS AND ACCESSORIES

Kindly refer to new basic seating benifit order form

### PLUS PACKAGES Select one below 15

**PLUS PACKAGE 1:**  
AADL CODE: W101  
CLIENT 250 - 349 LBS - \$900.00  
INCLUDES THE FOLLOWING:  
STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN  
SEAT WIDTH INCREASE FROM 16" TO 22"  
SEAT DEPTH INCREASE FROM 16" TO 24"  
5,6 OR 8 X 1 3/4" HD CASTERS  
HD STEEL FORKS

**PLUS PACKAGE 2:**  
AADL CODE: W102  
CLIENT 350 - 449 LBS - \$900.00  
INCLUDES THE FOLLOWING:  
STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN  
SEAT WIDTH INCREASE FROM 18" TO 26"  
SEAT DEPTH INCREASE FROM 16" TO 24"  
5,6 OR 8 X 1 3/4" HD CASTERS  
HD STEEL FORKS

**PLUS PACKAGE 3:**  
AADL CODE: W103  
CLIENT 450 - 600 LBS - \$1110.00  
INCLUDES THE FOLLOWING:  
STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN  
SEAT WIDTH INCREASE FROM 20" TO 30"  
SEAT DEPTH INCREASE FROM 16" TO 24"  
5,6 OR 8 X 1 3/4" HD CASTERS  
HD STEEL FORKS

70 Degree Pivot legrest are recommended for all PLUS packages. If 60 Degree or 80 Degree are required, Pin Style legrest are recommended.

### PLEASE NOTE:

- A1 =** Client's Cost
- A2 =** Clinical Justification
- A3 =** T-Arm recommended for seat depth & angle back adjustment.
- A4 =** To be mounted on inside of frame with inside wheels option.
- A5 =** Included when Heavy duty packages are ordered

**For more Seat To Floor combinations, visit our website**  
<https://powerplusbility.com/stf>  
**Please confirm with customer service ([csr@powerplusbility.com](mailto:csr@powerplusbility.com))**

### NOTES:

## INFORMATION

Customer Name: ..... Date: .....  
 Address: ..... PO#: .....  
 City: ..... Phone: .....  
 Account: .....  
 Authorizer Name: .....

BACKRESTS		O1
<input type="checkbox"/> N590	Relax Form Fit ( No Lateral )	\$690
<input type="checkbox"/> N590	Relax Gentle Fit ( <3" Lateral )	\$690
<input type="checkbox"/> N590	Relax Gentle Form Fit ( <3" Lateral )	\$690
<input type="checkbox"/> N590	Relax Gentle Spiral Fit ( <3" Lateral )	\$690
<input type="checkbox"/> N593	Relax Snug Form Fit ( >3" Lateral )	\$790
<input type="checkbox"/> N593	Relax Snug Spiral Fit ( >3" Lateral )	\$790
<input type="checkbox"/> N593	Relax Full Form Fit ( >3" Lateral )	\$790
<input type="checkbox"/> N593	Relax Full Spiral Fit ( >3" Lateral )	\$790
<input type="checkbox"/> N593	Relax Snug Fit ( >3" Lateral )	\$790
<input type="checkbox"/> N593	Relax Full Fit ( >3" Lateral )	\$790
Dimensions (in inches):		
Width: _____ x Height: _____		

ADDITIONAL OPTIONS		O2
<input type="checkbox"/> N500	Heavy-Duty Backrest ( >20" Wide )	\$300
<input type="checkbox"/> V571	Replacement Back Cover	\$183

PRESSURE REDUCTION OPTIONS		O3
Limited to 3 pieces per backrest.		
<input type="checkbox"/> V853	Lumbar Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/> V854	Lateral Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/> V855	Pelvic Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/> V856	Wedge Pressure Reduction QTY: ____ Location: ____	\$48.45
<input type="checkbox"/> V901	Thoracic Gel Pressure Reduction QTY: ____ Location: ____	\$69
<input type="checkbox"/> V902	Full Spinal Gel Pressure Reduction QTY: ____ Location: ____	\$103.50
<input type="checkbox"/> V903	Full Back Gel Pressure Reduction QTY: ____ Location: ____	\$137.08

CUSHIONS		O4
<input type="checkbox"/> V980	Relax Comfy ( Up to 20"x20" )	\$690
<input type="checkbox"/> V980HD	Relax Firm ( Over 20"x20" )	\$690
<input type="checkbox"/> V980	Relax AIRFLOW ( Up to 20"x20" )	\$790
<input type="checkbox"/> V980HD	Relax AIRFLOW ( Over 20"x20" )	\$790
Dimensions (in inches):		
Width: _____ x Depth: _____		

ADDITIONAL OPTIONS		O5
<input type="checkbox"/> V101	Replacement Cushion Cover	\$165
<input type="checkbox"/> V006	Incontinence Cover	\$135

HEADREST PAD		O7
<input type="checkbox"/> V081	Standard Pad	\$190
<input type="checkbox"/> N582	Headrest with Lateral Support	\$255

HEADREST HARDWARE		O8
<input type="checkbox"/> V338	Standard Straight Hardware	\$165
<input type="checkbox"/> N533	Multi-Axis Hardware ( L-Bend )	\$200
<input type="checkbox"/> N535	Gooseneck Multi-Axis ( S-Bend )	\$280
<input type="checkbox"/> V085	Mounting Plate <input type="checkbox"/> 1/2" or <input type="checkbox"/> 3/4"	\$189.70
<input type="checkbox"/> V570	Replacement Headrest Cover	\$120

ACCESSORIES		O6
<input type="checkbox"/> V577	Arm Trough Tray Assembly	\$432
<input type="checkbox"/> V476	Wheelchair Tray ( Up to 18" Wide )	\$257.60
<input type="checkbox"/> V542	Wheelchair Tray ( Over 18" Wide )	\$310.20
<input type="checkbox"/> V762	Elbow Stops	\$90
<input type="checkbox"/> V766	Lateral Support	\$90
<input type="checkbox"/> V769	Hand Block	\$90
<input type="checkbox"/> N814	Lateral Thigh Support	\$150
<input type="checkbox"/> V619	Modified Foot Support ( Padded Footbox )	\$359.10
<input type="checkbox"/> V581	Amputee Board ( Up to 20" Wide )	\$350
<input type="checkbox"/> V583	Heavy Duty Amputee Board	\$412.50
<input type="checkbox"/> V579	Pad for Amputee Board	\$62.50
<input type="checkbox"/> V572	Modular Arm Support ( Left Side )	\$62.50
<input type="checkbox"/> V572	Modular Arm Support ( Right Side )	\$62.50

**NOTES:**

