

STP STP Wheelchair (Tilt Wheelchair)



POWER PLUS MOBILITY

We will truly move you

Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____



FRAME COLOUR *Select one below* **1**

<input type="checkbox"/>	STPB	Tilt Frame - Black	2499
--------------------------	------	--------------------	------

CLIENT WEIGHT *Select one below* **2**

<input type="checkbox"/>	STPWA	Up to 250Lbs	N/C
<input type="checkbox"/>	STPWAMF	Exceeds 250Lbs	730
<input type="checkbox"/>	STPWAMG	Exceeds 350Lbs	900
<input type="checkbox"/>	STPWAMF	Exceeds 450Lbs	1110

SEAT WIDTH *Select one below* **3**

<input type="checkbox"/>	STPSW14	Seat Width 14"	N/C
<input type="checkbox"/>	STPSW15	Seat Width 15"	N/C
<input type="checkbox"/>	STPSW16	Seat Width 16"	N/C
<input type="checkbox"/>	STPSW17	Seat Width 17"	N/C
<input type="checkbox"/>	STPSW18	Seat Width 18"	N/C
<input type="checkbox"/>	STPSW19	Seat Width 19"	200
<input type="checkbox"/>	STPSW20	Seat Width 20"	200
<input type="checkbox"/>	STPSW21-30	Specify _____	450

21-30" Width

SEAT DEPTH *Select one below* **4**

<input type="checkbox"/>	STPSD16	Seat Depth 16"	N/C
<input type="checkbox"/>	STPSD17	Seat Depth 17"	N/C
<input type="checkbox"/>	STPSD18	Seat Depth 18"	N/C
<input type="checkbox"/>	STPSD19	Seat Depth 19"	200
<input type="checkbox"/>	STPSD20	Seat Depth 20"	200
<input type="checkbox"/>	STPSD21-24	Specify _____	450

21-24" Depth

FOOT REST AND FOOT PLATES *Select one below* **5**

<input type="checkbox"/>	STPOFST	Offset Foot Rest	N/C
<input type="checkbox"/>	STPF6A	60 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/>	STPF7A	70 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/>	STPF8A	80 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/>	STPEA	Elevating Legrest <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	404
<input type="checkbox"/>	STPAL	Articulating L/R <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	450
<input type="checkbox"/>	STPEB	Elevating Legrest System	450
<input type="checkbox"/>	STPPPL	Plastic Footplates	N/C
<input type="checkbox"/>	STPPXL	Extra Large Footplates	150
<input type="checkbox"/>	STPPAF	Angle Adjustable Footplates	150
<input type="checkbox"/>	STP1FB	One Piece Foot Board Padded	80

BACK TYPE **6**

<input type="checkbox"/>	STPAB1	Angle Adjustable Back	N/C
<input type="checkbox"/>	STPENDO	Endomorph Canes W ()	N/C
<input type="checkbox"/>	STPRBK	Recline Kit	560

Specify Width

ARM REST *Select one below* **7**

<input type="checkbox"/>	STPTA	T-Arm Ht. Adj./ Removable	N/C
<input type="checkbox"/>	STPFA	Flip Back Ht. Adj./ Removable	N/C
<input type="checkbox"/>	STPTAS	T-Arm Ht. Adj./ Removable	N/C
<input type="checkbox"/>	STPFAS	Flip Back Ht. Adj./ Removable	N/C
<input type="checkbox"/>	STPGWAPF	Gel Waterfall Arm Pads (A1)	150

Std 10"-14"

Std 10"-14"

Short 7"-11"

Short 7"-11"

STANDARD NO CHARGE SPECIFICATIONS

SEAT WIDTH: 14" – 18"
SEAT DEPTH (ADJUSTABLE)
SEAT HEIGHT (ADJUSTABLE): 10.25" – 18.25"
SEAT ANGLE: 0 – 30 [DEGREES]
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]
WEIGHT CAPACITY: 250 LBS
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)
TILT ADJUSTMENTS 0 TO 50 (DEGREES)
TWO-WAY LOCKING GAS STRUT
SINGLE RELEASE, DUAL CYLINDERS
STANDARD LOCKOUT COLLARS
REAR HEIGHT ADJUSTABLE ANTI TIPPERS

SEAT TO FLOOR SPECIFICATIONS

(LOW SEAT SETTING)			
STF	R.WHEEL	CASTER	FORK
<input type="checkbox"/>	S 10.50"	16"	3" XL SHORT
<input type="checkbox"/>	SL 13.00"	20"	5" XL SHORT
<input type="checkbox"/>	LL 13.75"	20"	5" XL SHORT
<input type="checkbox"/>	L 14.50"	20"	6" SHORT
<input type="checkbox"/>	L 14.50"	22"	5" SHORT
<input type="checkbox"/>	L14.50"	22"	6" SHORT
<input type="checkbox"/>	LH 15.25"	20"	5" SHT/LONG
<input type="checkbox"/>	LH 15.25"	20"	6" SHT/LONG
<input type="checkbox"/>	LH 15.25"	22"	5" SHT/LONG
<input type="checkbox"/>	LH 15.25"	22"	6" SHT/LONG
<input type="checkbox"/>	SH 16.00"	22"	6" LONG
(STANDARD SEAT SETTING)			
<input type="checkbox"/>	MH 16.75"	22"	6" LONG
<input type="checkbox"/>	H 17.50"	22"	6" LONG
<input type="checkbox"/>	H 17.50"	24"	6" LONG
<input type="checkbox"/>	A 18.25"	24"	8" LONG

Please Note : **A1** = Clients Responsibility

HEADREST		Select one below	8	PLUS PACKAGES		Select one below	15
<input type="checkbox"/>	STPMAHKP	Multi Axis Complete With Medium Pad	N/C	<input type="checkbox"/>	PLUS PACKAGE 1:		
<input type="checkbox"/>	STPMAHF	Multi Axis Complete With Large Pad	N/C		CLIENT 250 TO 349 LBS - \$730.00		
<input type="checkbox"/>	STPMAHFLAT	Multi Axis Complete Flat Pad (malleable)	100		INCLUDES THE FOLLOWING:		
REAR WHEELS				Select one below			
				9			
<input type="checkbox"/>	STPRWT0	20" Solid Treaded Tire	N/C		STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN		
<input type="checkbox"/>	STPRWT2	22" Solid Treaded Tire	N/C		SEAT WIDTH INCREASE FROM 16" TO 22"		
<input type="checkbox"/>	STPRWT4	24" Solid Treaded Tire	N/C		SEAT DEPTH INCREASE FROM 16" TO 24"		
STF 17-19"	<input type="checkbox"/>	STPRWT12I	12" Solid Tire Inside Mount	120	5,6 OR 8 X 1 3/4" HD CASTERS		
<input type="checkbox"/>	STPRWT16	16" Solid Tire	N/C	<input type="checkbox"/>	PLUS PACKAGE 2:		
<input type="checkbox"/>	STPAQR	Axles – Quick Release	120		CLIENT 350 TO 449 LBS - \$900.00		
<input type="checkbox"/>	STPHRT	Narrow Handrim Config. (1/2" less)	N/C		INCLUDES THE FOLLOWING:		
<input type="checkbox"/>	STPPCHR	Plastic Coated Hand Rims	N/C		STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN		
				SEAT WIDTH INCREASE FROM 18" TO 26"			
				SEAT DEPTH INCREASE FROM 16" TO 24"			
				5,6 OR 8 X 1 3/4" HD CASTERS			
				HD STEEL FORKS			
CASTER				Select one below			
				10			
<input type="checkbox"/>	STPFC1D	3" X 1" Caster	N/C	<input type="checkbox"/>	PLUS PACKAGE 3:		
<input type="checkbox"/>	STPFC1A	5" X 1.25" Caster	N/C		CLIENT 450 TO 600 LBS - \$1110.00		
<input type="checkbox"/>	STPFC1B	6" X 1.25" Caster	N/C		INCLUDES THE FOLLOWING:		
<input type="checkbox"/>	STPFC1C	8" X 1.25" Caster	N/C		STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN		
<input type="checkbox"/>	STPFC2A	5" X 1.75" Caster (Soft Roll)	80		SEAT WIDTH INCREASE FROM 20" TO 30"		
<input type="checkbox"/>	STPFC2B	6" X 1.75" Caster	80		SEAT DEPTH INCREASE FROM 16" TO 24"		
<input type="checkbox"/>	STPFC2C	8" X 1.75" Caster	80		5,6 OR 8 X 1 3/4" HD CASTERS		
<input type="checkbox"/>	STPFC4B	6" X 2 Caster (Steel Fork)	80		HD STEEL FORKS		
<input type="checkbox"/>	STPFC4C	8" X 2 Caster (Steel Fork)	80				
WHEEL LOCK				Select one below			
				11			
<input type="checkbox"/>	STPALPH	Aluminium Push To Lock	N/C	Notes:			
<input type="checkbox"/>	STPALPL	Aluminium Pull To Lock	N/C				
<input type="checkbox"/>	STPSTPH	Steel Push To Lock	N/C				
<input type="checkbox"/>	STPSTPL	Steel Pull To Lock	N/C				
<input type="checkbox"/>	STPBRXC	4" Brake Extension	40				
<input type="checkbox"/>	STPBRXA	6" Brake Extension	40				
<input type="checkbox"/>	STPBRXB	9" Brake Extension	130				
SEATING - BACKS							
				12			
				W () x H ()			
<input type="checkbox"/>	STPBK01	RELAX Gentle Fit Back	N/C				
SEATING - CUSHION				Select one below			
				13			
				W () x H ()			
14-20" wide	<input type="checkbox"/>	STPC01	RELAX Soft Cushion (Up to 250 lbs)	375			
14-20" wide	<input type="checkbox"/>	STPC02	RELAX Comfy Cushion (up to 250 lbs)	575			
21-28" wide	<input type="checkbox"/>	STPC03	RELAX Firm Cushion (up to 600 lbs)	575			
14-20" wide	<input type="checkbox"/>	STPC04	RELAX AIR FLOW Cushion	589			
ACCESSORIES AND POWER OPTION				14			
<input type="checkbox"/>	STPSBA	Seat Belt (Auto Style)	23				
<input type="checkbox"/>	STPCNC	Crutch And Cane Holder (A1)	60				
<input type="checkbox"/>	STPIVA	Adjustable I V Pole (A1)	130				
<input type="checkbox"/>	STPNLSP	Nylon Leg Strap Padded	50				
<input type="checkbox"/>	STPAMPL	Amputee Support Left	45				
<input type="checkbox"/>	STPAMPR	Amputee Support Right	45				
<input type="checkbox"/>	STPOXY	Oxygen Tank Holder	135				
<input type="checkbox"/>	STPTDFB	Transport Tie Downs Front	85				
<input type="checkbox"/>	STPTDRB	Transport Tie Downs Rear	85				
<input type="checkbox"/>	STPINTS	Independant Power Tilt System (A1)	2000				

Please Note : **A1** = Clients Responsibility

AADL CODE	PRODUCT	MSRP
SEATING OFF CHAIR		
BACKS		
N590	Relax Personal Fit - LESS THAN 3" LATERAL LEVEL B	\$690.00
N590	Relax Gentle Fit - LESS THAN 3" LATERAL LEVEL B	\$690.00
N593	Relax Snug Fit - GREATER THAN 3" LATERAL LEVEL B	\$790.00
N593	Relax Full Fit - GREATER THAN 3" LATERAL LEVEL B	\$790.00
N590	Relax Form Fit GENTLE - LESS THAN 3" LATERAL LEVEL B	\$690.00
N593	Relax Form Fit SNUG - GREATER THAN 3" LATERAL LEVEL B	\$790.00
N593	Relax Form Fit FULL - GREATER THAN 3" LATERAL LEVEL B	\$790.00
N590	Relax Spiral Fit GENTLE - LESS THAN 3" LATERAL LEVEL B	\$690.00
N593	Relax Spiral Fit SNUG - GREATER THAN 3" LATERAL LEVEL B	\$790.00
N593	Relax Spiral Fit FULL - GREATER THAN 3" LATERAL LEVEL B	\$790.00
N500	Oversized / Heavy Duty Option for Backrests - GREATER THAN 20" WIDE LEVEL B	\$300.00
V853	Lumbar/Extra Pressure Reduction - Maximum 3 pieces per back LEVEL A	\$48.45
V854	Lateral Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$48.45
V855	Pelvic Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$48.45
V856	Wedge Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$48.45
V901	Thoracic Gel Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$69.00
V902	Full Spinal Area Gel Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$103.50
V903	Full Back Area Gel Pressure Reduction - Maximum 3 Pieces Per Back LEVEL A	\$137.08
V571	REPLACEMENT BACK COVER FOR ALL SIZES LEVEL A OR LEVEL B	\$183.00
CUSHIONS		
V980	Relax Comfy CUSHION - Full Time User Sizes Up-to 20"x20" LEVEL A	\$589.00
V980 + V529	Relax Firm CUSHION - Full Time User Sizes Over 20"x20" (V529 = \$190) LEVEL A	\$779.00
V980	Relax AIRFLOW CUSHION - Adjustable Positioning Up-to 20"x20" LEVEL A	\$589.00
V980 + V529	Relax AIRFLOW CUSHION - Adjustable Positioning Over 20"x20" (V529 = \$190) LEVEL A	\$719.00
V101	REPLACEMENT CUSHION COVER FOR ALL SIZES LEVEL A OR LEVEL B	\$165.00
HEADRESTS		
V081	STANDARD PAD - F1 Contoured Pad Large LEVEL A	\$190.00
N582	HEADREST WITH LATERAL SUPPORT - FXL1 Flex Pad Large LEVEL B	\$255.00
V338	STANDARD STRAIGHT HARDWARE - LAH12 Vertical Horizontal Headrest Fixed 1/2" LEVEL A	\$165.00
N533	MULTI-AXIS HARDWARE - MAHL Multi-Axis Headrest (L-Bend) 1/2" LEVEL B	\$200.00
N535	GOOSENECK STYLE MULTI-AXIS - MAH Multi-Axis Headrest (S-Bend) 1/2" LEVEL B	\$280.00
V085	MOUNTING PLATE - UHRA Universal Headrest Plate - 1/2" OR 3/4" LEVEL A	\$189.70
V570	REPLACEMENT HEADREST COVER FOR ALL SIZES LEVEL A	\$120.00
ACCESSORIES		
V577	ARM TROUGH TRAY ASSEMBLY LEVEL A - CHPTL-P = Left Tray, CHPTR-P = Right Tray	\$432.00
V577	ARM TROUGH TRAY ASSEMBLY LEVEL A - Full padded tray with Gel Arm Through and Elbow Blocks	\$864.00
V476	WHEELCHAIR TRAY - Standard or Padded, Full or Half Up-to 18" Wide LEVEL A	\$257.60
V542	WHEELCHAIR TRAY - Standard or Padded, Full or Half Over 18" Wide LEVEL A	\$310.20
V762	ELBOW STOPS LEVEL A	\$90.00
V766	LATERAL SUPPORT LEVEL A	\$90.00
V769	HAND BLOCK LEVEL A	\$90.00
N814	LATERAL THIGH SUPPORT LEVEL B 150	\$150.00
V619-L	MODIFIED LEFT FOOT SUPPORT - Footbox, Padded Footbox, Contracture Sling (V619 \$179.55 X 2) LEVEL A	\$359.10
V619-R	MODIFIED RIGHT FOOT SUPPORT - Footbox, Padded Footbox, Contracture Sling (V619 \$179.55 X 2) LEVEL A	\$359.10
V581	AMPUTEE BOARD - Up-to 20" Wide LEVEL A	\$350.00
V583	HEAVY DUTY AMPUTEE BOARD - Up-to 28" x 22" LEVEL A	\$412.50
V579	PAD FOR AMPUTEE BOARD LEVEL A	\$62.50