

### HRZ Horizon Commode



**POWER PLUS MOBILITY**

*We will truly move you*

Customer Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ PO #: \_\_\_\_\_  
 City: \_\_\_\_\_ Phone #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Account #: \_\_\_\_\_ Authorizer Name: \_\_\_\_\_



**HORIZON 1** **1** OPTIONAL SEAT WIDTH **5**

- Weight Capacity upto 250 Lbs 1195
- 18" X 18" Urethane Seat
- Seat Height Including Seat Is 21"
- Seat width By Frame 18"
- Overall Width 20"
- Rear Axle Receiver Ready
- Stainless Steel Frame - Water Resistant
- 20" Back Height
- Black Urethane Easy To Clean Back Rest
- Black Molded Sanitary Cover
- Flip Back Arms
- Black Urethane Easy To Clean Arm Pads
- 6" X 10" X 14" Large Bucket
- 4" X 5" Rust Resistant Nylon Casters  
(3 Total Locking Casters & 1 Directional Lock Caster)

- HRZ14W 14" Seat Width 480
- HRZ15W 15" Seat Width 480
- HRZ16W 16" Seat Width 480
- HRZ17W 17" Seat Width 480
- HRZ23W 23" Seat Width 680
- HRZ24W 24" Seat Width 680
- HRZ25W 25" Seat Width 680
- HRZ26W 26" Seat Width 680

**HORIZON 2 - All options in HORIZON 1 plus:** **2**

- Weight Capacity upto 500 Lbs 1495
- 18" X 18" Easy to Clean Urethane Seat
- Seat Width By Frame 20"
- Overall With 22"
- Heavy Duty Construction

**HORIZON 3 - All options in HORIZON 1 plus:** **3**

- Weight Capacity upto 750 Lbs 1695
- 21" X 18" Urethane Seat
- Seat Width By Frame 22"
- Overall With 24"
- Extra Heavy Duty Construction

**OPTIONAL SEAT DEPTH** **4**


- HRZ14D 14" Seat Depth 480
- HRZ15D 15" Seat Depth 480
- HRZ16D 16" Seat Depth 480
- HRZ17D 17" Seat Depth 480
- HRZ18D 18" Seat Depth N/C
- HRZ19D 19" Seat Depth 480
- HRZ20D 20" Seat Depth 480
- HRZ21D 21" Seat Depth 680
- HRZ22D 22" Seat Depth 680
- HRZ23D 23" Seat Depth 680
- HRZ24D 24" Seat Depth 680
- HRZ25D 25" Seat Depth 680
- HRZ26D 26" Seat Depth 680

**ADDITIONAL OPTIONS** **6**

- HRZBH Custom Back Height 240
- HRZWC Capacity 750 To 1000 Lbs 1000
- HRZTC Tension Backrest 240
- HRZPB Padded Backrest 350
- HRZFS Form Fit Snug Water Proof Backrest & Hardware 790
- HRZLD Lock Down Flip Back Arms 240
- HRZCOS Closed Oval Soft Seat 18" To 22" 300
- HRZKHS Key Hole Soft Seat 18" To 22" 300
- HRZ60F 60° Footrest 240
- HRZ70F 70° Footrest N/C
- HRZ80F 80° Footrest 240
- HRZ90F 90° Footrest 240
- HRZEL Elevating Legrest 480  
(Not Water Proof)
- HRZAAF Angle Adjustable Footplates 240  
(Not Water Proof)
- HRZOPF One Piece Footboard 120
- HRZSB Seat Belt (Plastic) 26
- HRZAC Height Adjustable Casters 480
- HRZLC 4 Locking Casters (No Directional Lock) N/C
- HRZSTF Special Seat To Floor Height Specify: \_\_\_\_\_ 480
- HRZWH1 24" Wheels with 5" Casters 480
- HRZSAL Laterals 480

**Notes:**

**NOTES** **7**



A	Back Height:	<input style="width: 50px;" type="text"/>	G	Height:	<input style="width: 50px;" type="text"/>
B	Arm Height:	<input style="width: 50px;" type="text"/>	H	Shoulder Width:	<input style="width: 50px;" type="text"/>
C	Lower Leg Length:	<input style="width: 50px;" type="text"/>			
D	Seat Depth:	<input style="width: 50px;" type="text"/>		Age:	<input style="width: 50px;" type="text"/>

E Chest Width:

F Seat Width:

Weight: