

MP5

POWER WHEELCHAIR



ADP CODE: WPA3WM005

The MP5 wheelchair provides power and effortless maneuverability. This chair offers fast, easy and safe transportation. It is ideal for users looking for a bit more independence in getting around.



Single Drive PNG Joystick



Adjustable Rear Suspension



On Board Circuit Breaker



Neutral Levers To
Manually Maneuver
Power Chair



Height And Depth
Adjustable Seat Package



Standard Flip Up Foot board

CUSTOMIZE ANYTHING

Power Plus Mobility's advanced machinery specializes in manufacturing custom parts to customers' specifications from either a drawing or a sample. The parts we manufacture can be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

MPS

 **POWER PLUS MOBILITY**
Welcome home

STANDARD FEATURES

- Agile & Versatile Mid Wheel Drive Design
- Active Anti-Tippers & Suspension System
- Tight Turning Radius Of 24"
- 350 lbs Load Capacity
- Fully Programmable P & G VSI Controller
- Standard Flip Up Foot board
- Compact Adjustable Height
- Ergonomically Designed Captain Seat
- Rehab Frame Also Available
- Flip Back Arms

| SPECIFICATIONS | MP5 |
|-------------------------|--|
| OVERALL LENGTH: | 42" |
| OVERALL WIDTH: | 25" |
| WEIGHT CAPACITY: | 350 lbs. |
| SEATING: | 18"/20" CAPTAIN SEAT. REHAB FRAME 14" TO 26" |
| BATTERY SPECS.: | 12V 50A x 2 PCS |
| CONTROLLER TYPE: | P & G FULL PROGRAMMABLE |
| CHARGER TYPE: | 5 AMP, OFF BOARD 120/240 VOLT |
| MAX SPEED: | 8 km/h OR 5 MPH |
| WEIGHT WITH BATTERY: | 215 lbs. |
| WEIGHT WITHOUT BATTERY: | 150 lbs. |
| DRIVE WHEEL: | 12" x 4" SEMI PNEUMATIC |
| FRONT CASTER: | 8" x 2" SEMI PNEUMATIC |
| BATTERY RANGE: | 35 km (22 mi) |
| MOTOR TYPE: | 4600 RPM 450W 32:1 |
| SUSPENSION: | REAR CASTER |
| TURNING RADIUS: | 24" |

MP5
POWER WHEELCHAIR



Learn more about this product by scanning the QR Code

