

VOYAGER PLUS

FOLDING WHEELCHAIR



ADP CODE: WM0001063

AADL CODE: W153

The VOYAGER PLUS is built to withstand strenuous circumstances structured with heavy duty steel frame and cross braces. The VOYAGER PLUS is easily customizable and is available with bariatric weight capacities.



Vertical/Horizontal Rear Axle
Adjustability



Wrap Around Arms



Quick Release Arms



Dual Cross Braces



Angle Adjustable Back Canes 0°, 5°, 10° and 15°



Height Adjustable Back Canes 15"-24"

CUSTOMIZE ANYTHING

Power Plus Mobility's advanced machinery specializes in manufacturing custom parts to customers' specifications from either a drawing or a sample. The parts we manufacture can be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

VOYAGER PLUS



STANDARD FEATURES

- Chromoly Frame & Cross Braces
- Standard 18" Back Height
- 20", 22" or 24" Urethane Molded Rear Wheels
- New Light Weight 20" To 22" or 24" Spoked Rear Wheels
- Extruded Aluminum Forks
- Chrome Hand Rims
- 5", 6" or 8" Casters
- Black Nylon Upholstery
- Powder Coated Frame
- Flip Back Height Adjustable Removable Arms
- New Flip Back Arms (Allows Angle Back Adjustments)
- Detachable Swing In/Out Footrest with Heel Loops
- Precision Bearings
- Molded Plastic Footplates
- Seat Width 14" To 36"
- Black Padded Arm Pads
- Anti- Tip Rollers
- 6" Brake Extensions
- Double Cross Brace
- Overall Wheelchair Width: 27" (18" Seat Width Unfolded)
18" (18" Seat Width folded)

SPECIFICATIONS	VOYAGER PLUS
SEAT WIDTH:	14" TO 36"
SEAT DEPTH:	14" TO 24"
SEAT HEIGHT (ADJUSTABLE):	12.25" TO 19"
REAR WHEEL:	20", 22" OR 24" COMPOSITE URETHANE
CASTER SIZE:	5", 6" OR 8" URETHANE
BACK STYLE:	19" TO 22" LOW, 22" TO 25" HIGH SLIP ON BACK UPHOLSTERY
FRONT RIGGING:	SWING AWAY, DUAL SWING AWAY OR ELEVATING LEG RESTS
BACK ANGLE (ADJUSTABLE):	0°, 5°, 10° OR 15°
BACK HEIGHT RANGE	15" TO 24"
WEIGHT CAPACITY:	UP TO 600 lbs.
WEIGHT OF CHAIR:	36 lbs.

VOYAGER PLUS
FOLDING WHEELCHAIR



Learn more about this product by scanning the QR Code

