

EXT Extreme Wheelchair (Tilt Wheelchair)



Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____

FRAME COLOUR	Select one below	1
<input type="checkbox"/> EXTB	Dynamic Frame - Black	2299

CLIENT WEIGHT	Select one below	2
<input type="checkbox"/> EXTWA	Up to 250Lbs	N/C
<input type="checkbox"/> EXTWAMF	Exceeds 250Lbs	730
<input type="checkbox"/> EXTWAMG	Exceeds 350Lbs	992
<input type="checkbox"/> EXTWAMF	Exceeds 450Lbs	1110

SEAT WIDTH	Select one below	3
<input type="checkbox"/> EXTSW16	Seat Width 16"	N/C
<input type="checkbox"/> EXTSW17	Seat Width 17"	N/C
<input type="checkbox"/> EXTSW18	Seat Width 18"	N/C
<input type="checkbox"/> EXTSW19	Seat Width 19"	150
<input type="checkbox"/> EXTSW20	Seat Width 20"	150
<input type="checkbox"/> EXTSW21-30	Specify _____ (Must include Plus Package 1,2 or 3)	

SEAT DEPTH	Select one below	4
<input type="checkbox"/> EXTS16	Seat Depth 16"	N/C
<input type="checkbox"/> EXTS17	Seat Depth 17"	N/C
<input type="checkbox"/> EXTS18	Seat Depth 18"	N/C
<input type="checkbox"/> EXTS19	Seat Depth 19"	150
<input type="checkbox"/> EXTS20	Seat Depth 20"	150
<input type="checkbox"/> EXTS21-24	Specify _____ (Must include Plus Package 1,2 or 3)	

FOOT REST AND FOOT PLATES	Select one below	5
<input type="checkbox"/> EXTF6A	60 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> EXTF7A	70 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> EXTF8A	80 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> EXTEA	Elevating Legrest <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	404
<input type="checkbox"/> EXTAL	Articulating Leg Rests <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	450
<input type="checkbox"/> EXTEB	Elevating Legrest System	450
<input type="checkbox"/> EXTPPL	Plastic	N/C
<input type="checkbox"/> EXTPXL	Extra Large	150
<input type="checkbox"/> EXTPAF	Angle Adjustable Footplates	150
<input type="checkbox"/> EXTOFST	Offset Foot Rest Receivers	N/C
<input type="checkbox"/> EXT1FB	One Piece Foot Board	80

BACK TYPE		6
<input type="checkbox"/> EXTAB1	Angle Adjustable Back	N/C
<input type="checkbox"/> EXTENDO	Endomorph Canes W ()	N/C
<input type="checkbox"/> EXTRBK	Recline Kit	560

ARM REST	Select one below	7
<input type="checkbox"/> EXTITA	T-Arm Ht. Adj./ Removable (Std 10"-14")	N/C
<input type="checkbox"/> EXTITAS	T-Arm Ht. Adj./ Removable (Short 7"-11")	N/C
<input type="checkbox"/> EXTGWAPS	Gel Waterfall Arm Pads	150

STANDARD NO CHARGE SPECIFICATIONS
SEAT WIDTH: 14" – 18"
SEAT DEPTH (ADJUSTABLE): 16" – 18"
SEAT HEIGHT (ADJUSTABLE): 10.25" – 18.25"
SEAT ANGLE: 0 – 30 [DEGREES]
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]
WEIGHT CAPACITY: 250 LBS
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)
DYNAMIC TILT ADJUSTMENTS 0 TO 25 (DEGREES)
TWO-WAY LOCKING GAS STRUT
SINGLE RELEASE, DUAL CYLINDERS
STANDARD LOCKOUT COLLARS

SEAT TO FLOOR SPECIFICATIONS				
(LOW SEAT SETTING)				
STF	R. WHEEL	CASTER	FORK	
<input type="checkbox"/> S 10.50"	16"	3"	XL SHORT	
<input type="checkbox"/> SL 12.25	20"	5"	XL SHORT	
<input type="checkbox"/> SL 13.00"	20"	5"	XL SHORT	
<input type="checkbox"/> LL 13.75"	20"	5"	XL SHORT	
<input type="checkbox"/> L 14.50"	20"	6"	SHORT	
<input type="checkbox"/> L 14.50"	22"	5"	SHORT	
<input type="checkbox"/> L14.50"	22"	6"	SHORT	
<input type="checkbox"/> LH 15.25"	20"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	20"	6"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	6"	SHT/LONG	
<input type="checkbox"/> SH 16.00"	22"	6"	LONG	
(STANDARD SEAT SETTING)				
<input type="checkbox"/> MH 16.75"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	22"	8"	LONG	
<input type="checkbox"/> H 17.50"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	24"	8"	LONG	
<input type="checkbox"/> A 18.25"	24"	8"	LONG	

HEADREST		Select one below	8	SEATING - BACKS		Select one below	13
<input type="checkbox"/>	EXTMAHFLAT	Multi Axis Complete Flat Pad	N/C	W () x H ()			
<input type="checkbox"/>	EXTMAHF	Multi Axis Complete With Large Pad	N/C	<input type="checkbox"/>	EXTBK01	RELAX Gentle Fit Back	N/C
<input type="checkbox"/>	EXTMAHKP	Multi Axis Complete With Small Pad	N/C	PLUS PACKAGES		Select one below	14
REAR WHEELS		Select one below	9	<input type="checkbox"/> PLUS PACKAGE 1: CLIENT 250 TO 349 LBS - \$730.00 INCLUDES THE FOLLOWING: STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN SEAT WIDTH INCREASE FROM 16" TO 22" SEAT DEPTH INCREASE FROM 16" TO 24" 5,6 OR 8 X 1 1/4" HD CASTERS <input type="checkbox"/> PLUS PACKAGE 2: CLIENT 350 TO 449 LBS - \$992.00 INCLUDES THE FOLLOWING: STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN SEAT WIDTH INCREASE FROM 18" TO 26" SEAT DEPTH INCREASE FROM 16" TO 24" 5,6 OR 8 X 1 3/4" HD CASTERS HD STEEL FORKS <input type="checkbox"/> PLUS PACKAGE 3: CLIENT 450 TO 600 LBS - \$1110.00 INCLUDES THE FOLLOWING: STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN SEAT WIDTH INCREASE FROM 20" TO 30" SEAT DEPTH INCREASE FROM 16" TO 24" 5,6 OR 8 X 1 3/4" HD CASTERS HD STEEL FORKS			
<input type="checkbox"/>	EXTRWT0	20" Solid Treaded Tire	N/C				
<input type="checkbox"/>	EXTRWT2	22" Solid Treaded Tire	N/C				
<input type="checkbox"/>	EXTRWT4	24" Solid Treaded Tire	N/C				
<input type="checkbox"/>	EXTRWT16	16" Solid Tire	N/C				
<input type="checkbox"/>	EXTAQR	Axles – Quick Release	120				
<input type="checkbox"/>	EXTHRT	Narrow Handrim Config. (1/2" less)	N/C				
CASTER		Select one below	10				
<input type="checkbox"/>	EXTFC1D	3" X 1" Caster	N/C				
<input type="checkbox"/>	EXTFC1A	5" X 1.25" Caster	N/C				
<input type="checkbox"/>	EXTFC1B	6" X 1.25" Caster	N/C				
<input type="checkbox"/>	EXTFC1C	8" X 1.25" Caster	N/C				
<input type="checkbox"/>	EXTFC2A	5" X 1.75" Caster (Soft Roll)	80				
<input type="checkbox"/>	EXTFC2B	6" X 1.75" Caster	80				
<input type="checkbox"/>	EXTFC2C	8" X 1.75" Caster	80				
<input type="checkbox"/>	EXTFC4B	6" X 2 Caster	80				
<input type="checkbox"/>	EXTFC4C	8" X 2 Caster	80				
WHEEL LOCK		Select one below	11				
<input type="checkbox"/>	EXTALPH	Aluminium Push To Lock	N/C				
<input type="checkbox"/>	EXTALPL	Aluminium Pull To Lock	N/C				
<input type="checkbox"/>	EXTEXTH	Steel Push To Lock	N/C				
<input type="checkbox"/>	EXTEXTL	Steel Pull To Lock	N/C				
<input type="checkbox"/>	EXTBRXC	4" Brake Extension	N/C				
<input type="checkbox"/>	EXTBRXA	6" Brake Extension	N/C				
<input type="checkbox"/>	EXTBRXB	9" Brake Extension	130				
ACCESSORIES			12				
<input type="checkbox"/>	EXTSBA	Seat Belt (Auto Style)	26				
<input type="checkbox"/>	EXTCNC	Crutch And Cane Holder (A1)	60				
<input type="checkbox"/>	EXTIVA	Adjustable I V Pole (A1)	130				
<input type="checkbox"/>	EXTNLSP	Nylon Leg Strap Padded (A1)	50				

Please Note : A1 = Clients Responsibility

Notes: