

STP STP Wheelchair (Tilt Wheelchair)



Customer Name: _____
 Address: _____ PO #: _____
 City: _____ Phone #: _____ Date: _____
 Account #: _____ Authorizer Name: _____

FRAME COLOUR	Select one below	1
<input type="checkbox"/> STPB	Tilt Frame - Black	2299

CLIENT WEIGHT	Select one below	2
<input type="checkbox"/> STPWA	Up to 250Lbs	N/C
<input type="checkbox"/> STPWAMF	Exceeds 250Lbs	730
<input type="checkbox"/> STPWAMG	Exceeds 350Lbs	992
<input type="checkbox"/> STPWAMF	Exceeds 450Lbs	1110

SEAT WIDTH	Select one below	3
<input type="checkbox"/> STPSW16	Seat Width 16"	N/C
<input type="checkbox"/> STPSW17	Seat Width 17"	N/C
<input type="checkbox"/> STPSW18	Seat Width 18"	N/C
<input type="checkbox"/> STPSW19	Seat Width 19"	150
<input type="checkbox"/> STPSW20	Seat Width 20"	150
<input type="checkbox"/> STPSW21-30	Specify _____ (Must include Plus Package 1, 2 or 3)	

SEAT DEPTH	Select one below	4
<input type="checkbox"/> STPSD16	Seat Depth 16"	N/C
<input type="checkbox"/> STPSD17	Seat Depth 17"	N/C
<input type="checkbox"/> STPSD18	Seat Depth 18"	N/C
<input type="checkbox"/> STPSD19	Seat Depth 19"	150
<input type="checkbox"/> STPSD20	Seat Depth 20"	150
<input type="checkbox"/> STPSD21-24	Specify _____ (Must include Plus Package 1, 2 or 3)	

FOOT REST AND FOOT PLATES	Select one below	5
<input type="checkbox"/> STPOFST	Offset Foot Rest	N/C
<input type="checkbox"/> STPF6A	60 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> STPF7A	70 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> STPF8A	80 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> STPEA	Elevating Legrest <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	404
<input type="checkbox"/> STPAL	Articulating L/R <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	450
<input type="checkbox"/> STPEB	Elevating Legrest System	450
<input type="checkbox"/> STPPPL	Plastic	N/C
<input type="checkbox"/> STPPXL	Extra Large	150
<input type="checkbox"/> STPPAF	Angle Adjustable Footplates	150
<input type="checkbox"/> STP1FB	One Piece Foot Board	80

BACK TYPE		6
<input type="checkbox"/> STPAB1	Angle Adjustable Back	N/C
<input type="checkbox"/> STPENDO	Endomorph Canes W ()	N/C
<input type="checkbox"/> STPRBK	Recline Kit	560

ARM REST	Select one below	9
<input type="checkbox"/> STPTA	T-Arm Ht. Adj./ Removable (Std 10"-14")	N/C
<input type="checkbox"/> STPFA	Flip Back Ht. Adj./ Removable (Std 10"-14")	N/C
<input type="checkbox"/> STPTAS	T-Arm Ht. Adj./ Removable (Short 7"-11")	N/C
<input type="checkbox"/> STPFAS	Flip Back Ht. Adj./ Removable (Short 7"-11")	N/C
<input type="checkbox"/> STPGWAPF	Gel Waterfall Arm Pads	150

STANDARD NO CHARGE SPECIFICATIONS
SEAT WIDTH: 14" – 18"
SEAT DEPTH (ADJUSTABLE)
SEAT HEIGHT (ADJUSTABLE): 10.25" – 18.25"
SEAT ANGLE: 0 – 30 [DEGREES]
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]
WEIGHT CAPACITY: 250 LBS
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)
TILT ADJUSTMENTS 0 TO 50 (DEGREES)
TWO-WAY LOCKING GAS STRUT
SINGLE RELEASE, DUAL CYLINDERS
STANDARD LOCKOUT COLLARS

Tilt Options		8
<input type="checkbox"/> STPINTS	Independant Power Tilt System	2000

SEAT TO FLOOR SPECIFICATIONS				
(LOW SEAT SETTING)				
STF	R.WHEEL	CASTER	FORK	
<input type="checkbox"/> S 10.50"	16"	3"	XL SHORT	
<input type="checkbox"/> SL 12.25	20"	5"	XL SHORT	
<input type="checkbox"/> SL 13.00"	20"	5"	XL SHORT	
<input type="checkbox"/> LL 13.75"	20"	5"	XL SHORT	
<input type="checkbox"/> L 14.50"	20"	6"	SHORT	
<input type="checkbox"/> L 14.50"	22"	5"	SHORT	
<input type="checkbox"/> L14.50"	22"	6"	SHORT	
<input type="checkbox"/> LH 15.25"	20"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	20"	6"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	6"	SHT/LONG	
<input type="checkbox"/> SH 16.00"	22"	6"	LONG	
(STANDARD SEAT SETTING)				
<input type="checkbox"/> MH 16.75"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	22"	8"	LONG	
<input type="checkbox"/> H 17.50"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	24"	8"	LONG	
<input type="checkbox"/> A 18.25"	24"	8"	LONG	

HEADREST		Select one below	7	SEATING - BACKS		Select one below	14	
<input type="checkbox"/>	STPMAHFLAT	Multi Axis Complete Flat Pad	N/C	W () x H ()				
<input type="checkbox"/>	STPMAHF	Multi Axis Complete With Large Pad	N/C	<input type="checkbox"/>	STPBK01	RELAX Gentle Fit Back	N/C	
<input type="checkbox"/>	STPMAHKP	Multi Axis Complete With Small Pad	N/C	PLUS PACKAGES		Select one below	15	
REAR WHEELS		Select one below	10	<input type="checkbox"/> PLUS PACKAGE 1: CLIENT 250 TO 349 LBS - \$730.00 INCLUDES THE FOLLOWING: STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN SEAT WIDTH INCREASE FROM 16" TO 22" SEAT DEPTH INCREASE FROM 16" TO 24" 5,6 OR 8 X 1 3/4" HD CASTERS <input type="checkbox"/> PLUS PACKAGE 2: CLIENT 350 TO 449 LBS - \$992.00 INCLUDES THE FOLLOWING: STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN SEAT WIDTH INCREASE FROM 18" TO 26" SEAT DEPTH INCREASE FROM 16" TO 24" 5,6 OR 8 X 1 3/4" HD CASTERS HD STEEL FORKS <input type="checkbox"/> PLUS PACKAGE 3: CLIENT 450 TO 600 LBS - \$1110.00 INCLUDES THE FOLLOWING: STEEL FRAME WITH HEAVY DUTY BACK AND SEAT PAN SEAT WIDTH INCREASE FROM 20" TO 30" SEAT DEPTH INCREASE FROM 16" TO 24" 5,6 OR 8 X 1 3/4" HD CASTERS HD STEEL FORKS				
<input type="checkbox"/>	STPRWT0	20" Solid Treaded Tire	N/C					
<input type="checkbox"/>	STPRWT2	22" Solid Treaded Tire	N/C					
<input type="checkbox"/>	STPRWT4	24" Solid Treaded Tire	N/C					
(STF 17-19")	<input type="checkbox"/>	STPRWT12I	12" Solid Tire Inside Mount					120
<input type="checkbox"/>	STPRWT16	16" Solid Tire	N/C					
<input type="checkbox"/>	STPAQR	Axles – Quick Release	120					
<input type="checkbox"/>	STPHRT	Narrow Handrim Config. (1/2" less)	N/C					
CASTER		Select one below	11					
<input type="checkbox"/>	STPFC1D	3" X 1" Caster	N/C					
<input type="checkbox"/>	STPFC1A	5" X 1.25" Caster	N/C					
<input type="checkbox"/>	STPFC1B	6" X 1.25" Caster	N/C					
<input type="checkbox"/>	STPFC1C	8" X 1.25" Caster	N/C					
<input type="checkbox"/>	STPFC2A	5" X 1.75" Caster (Soft Roll)	80					
<input type="checkbox"/>	STPFC2B	6" X 1.75" Caster	80					
<input type="checkbox"/>	STPFC2C	8" X 1.75" Caster	80					
<input type="checkbox"/>	STPFC4B	6" X 2 Caster (Steel Fork)	80					
<input type="checkbox"/>	STPFC4C	8" X 2 Caster (Steel Fork)	80					
WHEEL LOCK		Select one below	12					
<input type="checkbox"/>	STPALPH	Aluminium Push To Lock	N/C					
<input type="checkbox"/>	STPALPL	Aluminium Pull To Lock	N/C					
<input type="checkbox"/>	STPSTPH	Steel Push To Lock	N/C					
<input type="checkbox"/>	STPSTPL	Steel Pull To Lock	N/C					
<input type="checkbox"/>	STPBRXC	4" Brake Extension	N/C					
<input type="checkbox"/>	STPBRXA	6" Brake Extension	N/C					
<input type="checkbox"/>	STPBRXB	9" Brake Extension	130					
ACCESSORIES			13					
<input type="checkbox"/>	STPSBA	Seat Belt (Auto Style)	26					
<input type="checkbox"/>	STPCNC	Crutch And Cane Holder (A1)	60					
<input type="checkbox"/>	STPIVA	Adjustable I V Pole (A1)	130					
<input type="checkbox"/>	STPNLSP	Nylon Leg Strap Padded (A1)	50					

Please Note : A1 = Clients Responsibility

Notes: