

STP STP Wheelchair (Tilt Wheelchair)



Customer Name: \_\_\_\_\_  
 Address: \_\_\_\_\_ PO #: \_\_\_\_\_  
 City: \_\_\_\_\_ Phone #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Account #: \_\_\_\_\_ Authorizer Name: \_\_\_\_\_



FRAME COLOUR	Select one below	1
<input type="checkbox"/> STPB	Tilt Frame - Black	2699
<input type="checkbox"/> STPS	Tilt Frame - Silver	2699
<input type="checkbox"/> STPC	Tilt Frame - Custom : _____	3099

CLIENT WEIGHT	Select one below	2
<input type="checkbox"/> STPWA	Up to 250lbs	N/C
<input type="checkbox"/> STPWAMF	Exceeds 250lbs	730
<input type="checkbox"/> STPWAMG	Exceeds 350lbs	992
<input type="checkbox"/> STPWAMH	Exceeds 400lbs (Max 600lbs)	1110

SEAT WIDTH	Select one below	3
<input type="checkbox"/> STPSW16	Seat Width 16"	N/C
<input type="checkbox"/> STPSW17	Seat Width 17"	N/C
<input type="checkbox"/> STPSW18	Seat Width 18"	N/C
<input type="checkbox"/> STPSW19	Seat Width 19"	399
<input type="checkbox"/> STPSW20	Seat Width 20"	399
<input type="checkbox"/> STPSWC	Custom Seat Width: _____	399

SEAT DEPTH	Select one below	4
<input type="checkbox"/> STPSD16	Seat Depth 16"	N/C
<input type="checkbox"/> STPSD17	Seat Depth 17"	N/C
<input type="checkbox"/> STPSD18	Seat Depth 18"	N/C
<input type="checkbox"/> STPSD19	Seat Depth 19"	265
<input type="checkbox"/> STPSD20	Seat Depth 20"	265
<input type="checkbox"/> STPSDC	Custom Seat Depth: _____	265

FOOT REST AND FOOT PLATES	Select one below	5
<input type="checkbox"/> STPOFST	Offset Foot Rest	N/C
<input type="checkbox"/> STPF6A	60 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> STPF7A	70 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> STPF8A	80 Degree Hanger <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	N/C
<input type="checkbox"/> STPEA	Elevating Legrest <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	404
<input type="checkbox"/> STPAL	Articulating L/R <input type="checkbox"/> Pivot <input type="checkbox"/> Pin	450
<input type="checkbox"/> STPEB	Elevating Legrest System	450
<input type="checkbox"/> STPPPL	Plastic	N/C
<input type="checkbox"/> STPPXL	Extra Large	150
<input type="checkbox"/> STPPAF	Angle Adjustable Footplates	240
<input type="checkbox"/> STPFRCH	Padded Footplates	88
<input type="checkbox"/> STP1FB	One Piece Foot Board	80
<input type="checkbox"/> STPNLS	Nylon Leg Strap Padded	60
<input type="checkbox"/> STPNLSC	Nylon Leg Strap Custom Size ( )	120
<input type="checkbox"/> STPCSL	Contracture Sling	240

STANDARD NO CHARGE SPECIFICATIONS
SEAT WIDTH: 14" – 18"
SEAT DEPTH (ADJUSTABLE)
SEAT HEIGHT (ADJUSTABLE): 12.25" – 18.25"
SEAT ANGLE: 0 – 30 [DEGREES]
BACK ANGLE (ADJUSTABLE): 0, 5, 10, 15 [DEGREES]
ADJUSTABLE STROLLER HANDLE: 0 – 180 [DEGREES]
WEIGHT CAPACITY: 250 LBS
WEIGHT OF CHAIR: 60 LBS (APPROXIMATELY)
TILT ADJUSTMENTS 0 TO 50 (DEGREES)
TWO-WAY LOCKING GAS STRUT
SINGLE RELEASE, DUAL CYLINDERS
STANDARD LOCKOUT COLLARS

Tilt Options	6	
<input type="checkbox"/> STPINTS	Independent Power Tilt System	2000

SEAT TO FLOOR SPECIFICATIONS				
(LOW SEAT SETTING)				
STF	R.WHEEL	CASTER	FORK	
<input type="checkbox"/> SL 12.25	20"	5"	XL SHORT	
<input type="checkbox"/> SL 13.00"	20"	5"	XL SHORT	
<input type="checkbox"/> LL 13.75"	20"	5"	XL SHORT	
<input type="checkbox"/> L 14.50"	20"	6"	SHORT	
<input type="checkbox"/> L 14.50"	22"	5"	SHORT	
<input type="checkbox"/> L14.50"	22"	6"	SHORT	
<input type="checkbox"/> LH 15.25"	20"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	20"	6"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	5"	SHT/LONG	
<input type="checkbox"/> LH 15.25"	22"	6"	SHT/LONG	
<input type="checkbox"/> SH 16.00"	22"	6"	LONG	
(STANDARD SEAT SETTING)				
<input type="checkbox"/> MH 16.75"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	22"	8"	LONG	
<input type="checkbox"/> H 17.50"	22"	6"	LONG	
<input type="checkbox"/> H 17.50"	24"	8"	LONG	
<input type="checkbox"/> A 18.25"	24"	8"	LONG	

BACK TYPE		7	SEATING - BACKS		Select one below	13		
<input type="checkbox"/>	STPAB1	Angle Adjustable Back	N/C	W ( ) x H ( )				
Specify Width	<input type="checkbox"/>	STPENDO	Endomorph Canes	W ( )	<input type="checkbox"/>	STPBK01	RELAX Gentle Fit Back	544
<input type="checkbox"/>	STPRBK	Recline Kit	560	<input type="checkbox"/>	STPBK02	RELAX Full Fit Back	608	
HEADREST		Select one below	8	<input type="checkbox"/>	STPBK03	RELAX Snug Fit Back	700	
<input type="checkbox"/>	STPMAHFLAT	Multi Axis Complete Flat Pad	460	<input type="checkbox"/>	STPBK04	RELAX Form Fit Back No Laterals	820	
<input type="checkbox"/>	STPMAHF	Multi Axis Complete With Large Pad	460	<input type="checkbox"/>	STPBK06	RELAX Form Fit Back Gentle Laterals	820	
<input type="checkbox"/>	STPMAHKP	Multi Axis Complete With Small Pad	460	<input type="checkbox"/>	STPBK07	RELAX Form Fit Back Full Laterals	820	
<b>NEW</b>	<input type="checkbox"/>	STPMAHFX	Multi Axis FLEX Complex Complete Large	525	<input type="checkbox"/>	STPBK05	RELAX Form Fit Back Snug Laterals	820
<input type="checkbox"/>	STPUPHT221	Contoured Pad	200	<input type="checkbox"/>	STPBK10	RELAX Spiral Fit Back Gentle Laterals	832	
(1/2" Tubing)	<input type="checkbox"/>	STPUPHT218	Std. vertical horizontal assembly (1/2")	256	<input type="checkbox"/>	STPBK09	RELAX Spiral Fit Back Snug Laterals	832
(3/4" Tubing)	<input type="checkbox"/>	STPUPHT217	Std. vertical horizontal assembly (3/4")	256	<input type="checkbox"/>	STPBK11	RELAX Spiral Fit Back Full Laterals	832
<input type="checkbox"/>	STPUHA	Universal Hardware Adaptor	100	<input type="checkbox"/>	STPTAB	RELAX Tension Back Upto 18"	140	
<b>NEW</b>	<input type="checkbox"/>	STPMAHGA-P	Multi Axis with Headrest Plate	224	<input type="checkbox"/>	STPTAB+	RELAX Tension Back Over 18"	240
ARM REST		Select one below	9	<input type="checkbox"/>	STPBGL1	Half Overlay Gel (3 Piece)	348	
(Std 10"-14")	<input type="checkbox"/>	STPTA	T-Arm Ht. Adj./ Removable	240	<input type="checkbox"/>	STPBGL2	Full Overlay Gel (6 Pieces)	696
(Std 10"-14")	<input type="checkbox"/>	STPF A	Flip Back Ht. Adj./ Removable	340	<input type="checkbox"/>	STPPSTP	Pin Extension (1" Extension)	40
(Short 7"-11")	<input type="checkbox"/>	STPTAS	T-Arm Ht. Adj./ Removable	240	<input type="checkbox"/>	STPBK12	*Custom Width : _____	176
(Short 7"-11")	<input type="checkbox"/>	STPFAS	Flip Back Ht. Adj./ Removable	340	<input type="checkbox"/>	STPBK13	*Custom Height : _____	176
<input type="checkbox"/>	STPGWAPF	Gel Waterfall Arm Pads	180	<b>*Note: Please refer to the Seating Booklet for custom sizes</b>				
REAR WHEELS		Select one below	10	SEATING - CUSHION		Select one below	14	
<input type="checkbox"/>	STPRWT0	20" Solid Treaded Tire	N/C	W ( ) x D ( )				
<input type="checkbox"/>	STPRWT2	22" Solid Treaded Tire	N/C	<input type="checkbox"/>	STPC01	RELAX Soft Cushion	480	
<input type="checkbox"/>	STPRWT4	24" Solid Treaded Tire	N/C	<input type="checkbox"/>	STPC02	RELAX Comfy Cushion	624	
(STF 13-19")	<input type="checkbox"/>	STPRWT12	12" Solid Tire	N/C	<input type="checkbox"/>	STPC03	RELAX Firm Cushion	784
(STF 17-19")	<input type="checkbox"/>	STPRWT12I	12" Solid Tire Inside Mount	120	<input type="checkbox"/>	STPCGL1	Half Overlay Gel (1 Piece)	116
<input type="checkbox"/>	STPRWT16	16" Solid Tire	N/C	<input type="checkbox"/>	STPCGL2	Full Overlay Gel (2 Pieces)	232	
<input type="checkbox"/>	STPRWA2	22" Air Tire	120	<input type="checkbox"/>	STPCCL	Contoured Leg Channel (Bolster) Std.	180	
<input type="checkbox"/>	STPRWA4	24" Air Tire	120	<input type="checkbox"/>	STPCCLC	Contoured Leg Channel (Bolster) Custom	320	
<input type="checkbox"/>	STPAQR	Axles – Quick Release	120	<b>*Note: Please refer to the Seating Booklet for custom sizes</b>				
<input type="checkbox"/>	STPHRT	Narrow Handrim Config. (1/2" less)	N/C	<input type="checkbox"/>	STPC00	Extra Cushion Cover	99	
CASTER		Select one below	11	ACCESSORIES		15		
<input type="checkbox"/>	STPFC1A	5" X 1.25" Caster	N/C	<input type="checkbox"/>	STPSBA	Seat Belt (Auto Style)	26	
<input type="checkbox"/>	STPFC1B	6" X 1.25" Caster	N/C	<input type="checkbox"/>	STPSBAP	Seat Belt (Padded)	56	
<input type="checkbox"/>	STPFC1C	8" X 1.25" Caster	N/C	<input type="checkbox"/>	STPSBA4P	Seat Belt (4 Point Padded)	130	
<input type="checkbox"/>	STPFC2A	5" X 1.75" Caster (Soft Roll)	80	<input type="checkbox"/>	STPTR A	Full Padded Tray Locking	280	
<input type="checkbox"/>	STPFC2B	6" X 1.75" Caster	80	<input type="checkbox"/>	STPTRB	Full Padded Tray Slide On	280	
<input type="checkbox"/>	STPFC2C	8" X 1.75" Caster	80	<input type="checkbox"/>	STPTRC	Full Plastic Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	280	
<input type="checkbox"/>	STPFC4B	6" X 2 Caster (Steel Fork)	80	<input type="checkbox"/>	STPTRD	Half Padded Left Tray	240	
<input type="checkbox"/>	STPFC4C	8" X 2 Caster (Steel Fork)	80	<input type="checkbox"/>	STPTRE	Half Padded Right Tray	240	
WHEEL LOCK		Select one below	12	<input type="checkbox"/>	STPTRF	Half Left Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	240	
<input type="checkbox"/>	STPALPH	Aluminium Push To Lock	N/C	<input type="checkbox"/>	STPTRG	Half Right Tray <input type="checkbox"/> ABS Black <input type="checkbox"/> Clear	240	
<input type="checkbox"/>	STPALPL	Aluminium Pull To Lock	N/C	<input type="checkbox"/>	STPTDFB	Transport Tie Downs Front	85	
<input type="checkbox"/>	STPSTPH	Steel Push To Lock	N/C	<input type="checkbox"/>	STPTDRB	Transport Tie Downs Rear	85	
<input type="checkbox"/>	STPSTPL	Steel Pull To Lock	N/C	<input type="checkbox"/>	STPOXY	Oxygen Tank Holder	195	
<input type="checkbox"/>	STPBRXC	4" Brake Extension	N/C	<input type="checkbox"/>	STPIVA	Adjustable I V Pole	130	
<input type="checkbox"/>	STPBRXA	6" Brake Extension	N/C	<input type="checkbox"/>	STPCHPT-L	Left Arm Trough Tray	400	
<input type="checkbox"/>	STPBRXB	9" Brake Extension	130	<input type="checkbox"/>	STPCHPT-R	Right Arm Trough Tray	400	
<b>NEW</b>	<input type="checkbox"/>	STPATTBRX	Attendant Brakes	240	<input type="checkbox"/>	STPLATB	Lateral Bite out	180
				<input type="checkbox"/>	STPADD8	8" Adductor	164	
				<input type="checkbox"/>	STPADD10	10" Adductor	164	
				<input type="checkbox"/>	STPADD12	12" Adductor	164	