

# WHEELCHAIR BOOKLET

TILT WHEELCHAIRS  
STP  
EXTREME  
NV TILT

FOLDING WHEELCHAIRS  
MAGIC PLUS  
VOYAGER PLUS  
GLYDER

POWER WHEELCHAIR  
MP5

COMMODOES  
HORIZON NON-TILT  
HORIZON TILT

SIT BACK AND RELAX

# TABLE OF CONTENTS

## TILT WHEELCHAIR

01



STP

05



EXTREME

09



RV TILT

## FOLDING WHEELCHAIR

13



MAGIC PLUS

17



VOYAGER PLUS

21



GLYDER

## POWER WHEELCHAIR

25



MP5

## COMMODES

29



HORIZON

31



HORIZON TILT

*At Power Plus Mobility, We care!*

Thank you for considering a wheelchair from Power Plus Mobility Inc. We are always improving in our ways to serve you better and would love to hear your questions or comments about our literatures, the reliability of our chairs and service you received while making the decision of choosing our wheelchair. Please feel free to contact us via our email address: [customer@powerplusbility.com](mailto:customer@powerplusbility.com)

# STP

SUPER TILT PLUS

TILT WHEELCHAIR



One of Power Plus Mobility's best selling wheelchairs. The **STP** is a light weight tilt chair designed for adjustment and maneuverability. The **STP** tilt system provides a change of orientation and thus redistributes pressure from one area to another area.




UP TO 55° DYNAMIC TILT




Power Tilt Option



Center Pivot Swing In And Swing Out Front Rigging



Dual Cylinders Allows For Consistent Tilt For All Weight Capacities



Narrow Configuration  
(Inside Wheel Mount)



Rear Axle Adjustability



Rear Depth Adjustable Frame

## ABOUT STP

One of Power Plus Mobility's best selling wheelchairs. The **STP** is designed for adjustment and maneuverability. With comfort being another major point in this design, our **STP** will be your choice in tilt wheelchairs. We promise this product will exceed your expectations. By meeting all requirements and necessities for a wheelchair like this and even exceeding them. The **STP** offers a wide range of dynamics that allow this wheelchair to accommodate many different situations, assuring a right fit for all, even in the most difficult positions that other wheelchairs cannot serve.

## CUSTOMIZE ANYTHING

Power Plus Mobility advanced machinery specializes in manufacturing custom parts to customers specifications from either a drawing or a sample. The parts we manufacture could be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

**STP**  
SUPER TILT PLUS



## STANDARD FEATURES

- Light Weight Frame
- Depth Adjustable Frame
- ABS Seat Pan
- Dual Gas Struts For Easy Tilt (Up To 55° )
- Tilt Collar lock, locked out at 33°
- New Recessed Lower Frame For Foot Propulsion
- Dual Gas Cylinders for Smooth & Easy Tilting
- Soft Set Gas Cylinder Release Mechanism
- Silver Grey or Black Frame Colour
- T-Style Arms
- Adjustable Stroller Handle
- Swing In & Out Footrest
- Adjustable Caster Journals
- Recessed Axle Plates

SPECIFICATIONS	STP
SEAT WIDTH:	14" TO 36"
SEAT DEPTH:	16" TO 24"
SEAT HEIGHT:	(ADJUSTABLE): 12.25" TO 18.25"
TILT ANGLE:	0° TO 55°
BACK ANGLE:	(ADJUSTABLE): 0° , 5° , 10° OR 15°
ADJUSTABLE STROLLER HANDLE:	0° TO 180°
WEIGHT CAPACITY:	UP TO 600 lbs.
WEIGHT OF CHAIR:	60 lbs. (APPROXIMATELY)
	TWO-WAY LOCKING GAS STRUT
	SINGLE RELEASE, DUAL CYLINDERS
	STANDARD LOCKOUT COLLARS (33° )



PP-MWCSTP

# EXTREME

## TILT WHEELCHAIR



The **EXTREME** tilt Wheelchair from Power Plus Mobility has a light frame and is depth adjustable. The **EXTREME** tilt system provides a change of orientation and thus redistributes pressure from one area to another. The front seat to floor height also remains constant during 0 To 25 tilt degrees



Seat To Floor Height  
Adjustability




Zero Maintenance Soft Set Cylinders



Offset Front Hangers  
(Optional)





Dual Cylinders Allows For Consistent Tilt For All Weight Capacities



Standard Angle Adjustable Stroller Handles



Angle Adjustable Front Journals

## ABOUT EXTREME

Our wheelchair **EXTREME** lives up to its name. Exceeding the boundaries of what you would find in a wheelchair, the **EXTREME** is in a category all by itself. When Power Plus Mobility first began the process of creating **EXTREME**, it was definite that this product would deliver an extraordinary result to the customer. With its ability to be fully adjustable and easily customizable due to our new technology, the **EXTREME** is perfect for many. Designed to meet all the users needs, this wheelchair rises above others sent out by competitors. We built it with our customers in mind and every detail was fully analyzed and researched during the design process to ensure the most effective characteristics were chosen. The **EXTREME** is a prime example of the hard work we put into our products.

## CUSTOMIZE ANYTHING

Power Plus Mobility advanced machinery specializes in manufacturing custom parts to customers specifications from either a drawing or a sample. The parts we manufacture could be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

# EXTREME




## STANDARD FEATURES

- Light Weight Frame
- Depth Adjustable Frame
- ABS Seat Pan
- Dual Gas Struts for Easy Tilt (Up To 25° )
- New Recessed Lower Frame For Foot Propulsion
- Dual Gas Cylinders for Smooth & Easy Tilting
- Relax Seating System
- New Gas Cylinder Release Mechanism
- Black Frame Colour
- T-Style Arms
- Adjustable Stroller Handle
- Swing In & Out Footrest
- Adjustable Caster Journals
- New Recessed Axles
- Maintains Front STF Height In Tilt

SPECIFICATIONS	Extreme
SEAT WIDTH:	14" TO 36"
SEAT DEPTH:	(ADJUSTABLE): 16" TO 24"
SEAT HEIGHT:	(ADJUSTABLE): 12.25" TO 18.25"
Tilt ANGLE:	0° TO 25°
BACK ANGLE:	(ADJUSTABLE): 0°, 5°, 10° OR 15°
ADJUSTABLE STROLLER HANDLE:	0° TO 180°
WEIGHT CAPACITY:	UP TO 600 lbs.
WEIGHT OF CHAIR:	60 lbs. (APPROXIMATELY)
	DYNAMIC TILT ADJUSTMENTS 0° TO 25°
	TWO-WAY LOCKING GAS STRUT
	SINGLE RELEASE, DUAL CYLINDERS
	STANDARD LOCKOUT COLLARS

**EXTREME**

TILT WHEELCHAIR



PP-MWCEXT

# NV TILT


## TILT WHEELCHAIR




The **NV TILT** is a light weight tilt chair designed for adjustment and maneuverability. The **NV TILT** tilt system provides a change of orientation and thus redistributes pressure from one area to another area.



Center Pivot Swing In And Swing Out  
Front Rigging  
(Optional)



Dual Cylinders Allows For Consistent  
Tilt For All Weight Capacities



Narrow Configuration  
(Inside Wheel Mount)



Rear Axle Adjustability



Rear Depth Adjustable Frame

## ABOUT NV TILT

One of Power Plus Mobility's best selling wheelchairs. The **NV TILT** is designed for adjustment and maneuverability. With comfort being another major point in this design, our **NV TILT** will be your choice in tilt wheelchairs. We promise this product will exceed your expectations. By meeting all requirements and necessities for a wheelchair like this and even exceeding them. The **NV TILT** offers a wide range of dynamics that allow this wheelchair to accommodate many different situations, assuring a right fit for many, even in the most difficult positions that other wheelchairs cannot serve.

## CUSTOMIZE ANYTHING

Power Plus Mobility advanced machinery specializes in manufacturing custom parts to customers specifications from either a drawing or a sample. The parts we manufacture could be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

# NV TILT



## STANDARD FEATURES

- Light Weight Frame
- Depth Adjustable Frame
- ABS Seat Pan
- Dual Gas Struts For Easy Tilt (Up To 55°)
- Tilt Locked out at 33°
- New Recessed Lower Frame For Foot Propulsion
- Dual Gas Cylinders for Smooth & Easy Tilting
- Soft Set Gas Cylinder Release Mechanism
- Silver Grey or Black Frame Colour
- T-Style Arms
- Adjustable Stroller Handle
- Adjustable Caster Journals
- Recessed Axle Plates
- One Piece Elevating System
- Incontinent Seating

SPECIFICATIONS	STP
SEAT WIDTH:	14" TO 36"
SEAT DEPTH:	16" TO 24"
SEAT HEIGHT:	(ADJUSTABLE): 12.25" TO 18.25"
TILT ANGLE:	0° TO 55°
BACK ANGLE:	(ADJUSTABLE): 0° , 5° , 10° OR 15°
ADJUSTABLE STROLLER HANDLE:	0° TO 180°
WEIGHT CAPACITY:	UP TO 600 lbs.
WEIGHT OF CHAIR:	60 lbs. (APPROXIMATELY)
	TWO-WAY LOCKING GAS STRUT
	SINGLE RELEASE, DUAL CYLINDERS
	STANDARD LOCKOUT COLLARS (33°)




PP-MWCSTP

# MAGIC PLUS

FOLDING WHEELCHAIR



The **MAGIC PLUS** is a lightweight aluminum folding chair with vertical multi-axle positioning. This chair also offers one of the lowest seat to floor heights.



Vertical Rear Axle Adjustability  
Ranging From 12.5" To 18.25"



Angle Adjustable Back Canes 0°, 5°,  
10° Or 15°




Dual cross braces





Two Point Flip Back Arm



Height Adjustable Arm



Swing In And Swing Out Front Rigging

## ABOUT MAGIC PLUS

Power Plus Mobility's newly ADP certified **MAGIC PLUS** is one of our newest manual wheelchairs. Designed with a touch of magic and an aluminum lightweight composition, this accessible model will be useful for customers for everyday use. Here at Power Plus Mobility, we promise a product that will exceed your expectations with all requirements and necessities for a wheelchair like this and that will always be our goal. Add a little magic to your life and order today.

## CUSTOMIZE ANYTHING

Power Plus Mobility advanced machinery specializes in manufacturing custom parts to customers specifications from either a drawing or a sample. The parts we manufacture could be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

# MAGIC PLUS



## KEY FEATURES

- Light Weight Aluminum Frame & Cross Braces
- Multi Axle Position - Vertical
- 20", 22" or 24" Nylon Molded Rear Wheels
- New Light Weight 20", 22" or 24" Spoked Rear Wheels
- Extruded Aluminum Forks
- Chrome Hand Rims
- 5", 6" or 8" Casters
- Black Nylon Upholstery
- Powder Coated Frame
- T-Arm Height Adjustable/ Removable
- New Flip Back Arms Which Allow Angle Back Adjustment
- Detachable Swing In/Out Footrest with Heel Loops
- Precision Bearings
- Molded Plastic Footplates
- Seat Width 14" To 20"
- Black Padded Arm Pads
- Anti-Tip Rollers
- 6" Brake Extensions

SPECIFICATIONS	MAGIC PLUS
SEAT WIDTH:	14" TO 20"
SEAT DEPTH:	14" TO 24"
SEAT HEIGHT (ADJUSTABLE):	12.5" TO 18.25"
REAR WHEEL:	20",22" OR 24" COMPOSITE URETHANE
CASTER SIZE:	5", 6" OR 8" URETHANE
BACK STYLE:	19" TO 22" LOW, 22" TO 25" HIGH SLIP ON BACK UPHOLSTERY
FRONT RIGGING:	SWINGAWAY, DUAL SWINGAWAY OR ELEVATING LEGRESTS
BACK ANGLE (ADJUSTABLE):	0° , 5° , 10° OR 15°
BACK HEIGHT RANGE:	15" TO 24"
WEIGHT CAPACITY:	250 lbs.
WEIGHT OF CHAIR:	31 lbs.

MAGIC PLUS  
FOLDING WHEELCHAIR



PP-MWCMGP

# VOYAGER PLUS

FOLDING WHEELCHAIR



**VOYAGER PLUS** by Power Plus Mobility was built in mind to withstand strenuous circumstances. Structured with a lightweight chromoly frame and cross braces, it can be challenged with the reality of anyone's busy life in a reliable form. **VOYAGER PLUS** comes with bariatric weight capacities.



Vertical/Horizontal Rear Axle  
Adjustability



Wrap Around Arms



Quick Release Arms



Dual Cross Braces



Angle Adjustable Back Canes 0°, 5°, 10° and 15°



Height Adjustable Back Canes 15"-24"

## ABOUT VOYAGER PLUS

ADP approved, the **VOYAGER PLUS** by Power Plus Mobility was built in mind to withstand strenuous circumstances. Structured with a lightweight chromoly frame and cross braces, it can be challenged with the reality of anyone's busy life in a reliable form. This wheelchair is easily customizable to fit the needs required in order to bring comfortability, yet sturdiness and quality at the same time. **VOYAGER PLUS** is just one example of how our product's will continually and truly move you. If interested the information provided below can give you further insight into this certain product. You can guarantee that this product is made with the intention to deliver the maximum efficiency possible as are all Power Plus Mobility products.

## CUSTOMIZE ANYTHING

Power Plus Mobility advanced machinery specializes in manufacturing custom parts to customers specifications from either a drawing or a sample. The parts we manufacture could be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

# VOYAGER PLUS



## FEATURES

- Chromoly Frame & Cross Braces
- Standard 18" Back Height
- 20", 22" or 24" Nylon Molded Rear Wheels
- New Light Weight 20" To 22" or 24" Spoked Rear Wheels
- Extruded Aluminum Forks
- Chrome Hand Rims
- 5", 6" or 8" Casters
- Black Nylon Upholstery
- Powder Coated Frame
- Flip Back Height Adjustable Removable Arms
- New Flip Back Arms Which Allow Angle Back Adjustments
- Detachable Swing In/Out Footrest with Heel Loops
- Precision Bearings
- Molded Plastic Footplates
- Seat Width 14" To 36"
- Black Padded Arm Pads
- Anti- Tip Rollers
- 6" Brake Extensions
- Double Cross Brace

SPECIFICATIONS	VOYAGER PLUS
SEAT WIDTH:	14" TO 36"
SEAT DEPTH:	14" TO 24"
SEAT HEIGHT (ADJUSTABLE):	12.25" TO 19"
REAR WHEEL:	20", 22" OR 24" COMPOSITE URETHANE
CASTER SIZE:	5", 6" OR 8" URETHANE
BACK STYLE:	19" TO 22" LOW, 22" TO 25" HIGH SLIP ON BACK UPHOLSTERY
FRONT RIGGING:	SWINGAWAY, DUAL SWINGAWAY OR ELEVATING LEGRESTS
BACK ANGLE (ADJUSTABLE):	0°, 5°, 10° OR 15°
BACK HEIGHT RANGE	15 TO 24
WEIGHT CAPACITY:	UP TO 600 lbs.
WEIGHT OF CHAIR:	36 lbs.

**VOYAGER PLUS**  
FOLDING WHEELCHAIR




PP-MWCVOY

# GLYDER

FOLDING WHEELCHAIR



The **GLYDER** is structured with a lightweight aluminum frame and cross brace. Horizontally and vertically versatile, **GLYDER** is one of our newest wheelchairs, designed and constructed with accessibility and maximum mobility in any situation. Possessing a variety of customizable options.



Horizontal Rear Axle Adjustability  
Ranging From 12.5" To 18.25"



Angle Adjustable Back Canes 0°, 5°,  
10° Or 15°




Dual Cross Braces





Two Point Flip Back Arm



Height Adjustable Arm



Swing In And Swing Out Front Rigging

## ABOUT GLYDER

Power Plus Mobility's **GLYDER** is structured with a lightweight aluminum frame and cross brace. Horizontally and vertically versatile, **GLYDER** is one of our newest wheelchairs, designed and constructed with accessibility and maximum mobility in any situation. Possessing a variety of customizable options, **GLYDER** is capable of any situation for any individual. If interested the information provided below can give you further insight into this certain product. We promise that this product is made with the intention to deliver the maximum efficiency as possible, thus are all Power Plus Mobility products, where our products will truly move you.

## CUSTOMIZE ANYTHING

Power Plus Mobility advanced machinery specializes in manufacturing custom parts to customers specifications from either a drawing or a sample. The parts we manufacture could be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

# GLYDER




## KEY FEATURES

- Light Weight Aluminum Frame & Cross Braces
- Multi Axle Position - Vertical & Horizontal
- 20", 22" Or 24" Nylon Molded Rear Wheels
- New Light Weight 20", 22"Or 24" Spoked Rear Wheels
- Extruded Aluminum Forks
- Chrome Hand Rims
- 5", 6" Or 8" Casters
- Black Nylon Upholstery
- Powder Coated Frame
- T-Arm Height Adjustable/ Removable
- New Flip Back Arms Which Allow Angle Back Adjustments
- Detachable Swing In/Out Footrest with Heel Loops
- Precision Bearings
- Molded Plastic Footplates
- Seat Width 14" To 20" with Aluminum Frame & 21" > with Robust Steel Frame
- Black Padded Arm Pads
- Anti-Tip Rollers
- 6" Brake Extensions

SPECIFICATIONS	GLYDER
SEAT WIDTH:	14" TO 20" ALUMINUM FRAME 21"> ROBUST STEEL FRAME
SEAT DEPTH:	14" TO 24"
SEAT HEIGHT (ADJUSTABLE):	12.5" TO 18.25"
REAR WHEEL:	20",22" OR 24" COMPOSITE URETHANE
CASTER SIZE:	4",5",6" OR 8" URETHANE
BACK STYLE:	19" TO 22" LOW, 22" TO 25" HIGH SLIP ON BACK UPHOLSTERY
FRONT RIGGING:	SWINGAWAY, DUAL SWINGAWAY OR ELEVATING LEGRESTS
BACK ANGLE (ADJUSTABLE):	0° , 5° , 10° OR 15°
BACK HEIGHT RANGE:	15" TO 24"
WEIGHT CAPACITY:	250 lbs.
WEIGHT OF CHAIR:	31 lbs.

**GLYDER**  
FOLDING WHEELCHAIR



PP-MWCGLY

# MPS

POWER WHEELCHAIR



The **MPS** wheelchair will provide the power and maneuverability wherever it is used. It is an ideal product for those in need of a wheelchair that patients and/or customers can rely on to provide fast, easy, and especially safe transportation.




Single Drive PNG Joystick



Adjustable Rear Suspension



On Board Circuit Breaker



Neutral Levers To Manually Maneuver Power Chair



Height And Depth Adjustable Seat Package



Standard Flip Up Footboard

## ABOUT MP5

The Power wheelchair **MP5** is essentially a product from Power Plus Mobility that will surely deliver the power that it intends; after all it is in our name and we can fully promise that this product is the best it can be for our customers. This wheelchair will provide the power and maneuverability wherever it is used. It is an ideal product for those in need of a wheelchair that patients and/or customers can rely on to provide fast, easy, and especially safe transportation. If interested the information provided below can give you further insight into this certain product. You can guarantee that this product is made with intention to deliver the maximum efficiency as possible as are all Power Plus Mobility products.

## CUSTOMIZE ANYTHING

Power Plus Mobility advanced machinery specializes in manufacturing custom parts to customers specifications from either a drawing or a sample. The parts we manufacture could be simple or more complex and can be manufactured in small or large quantities. Contact your local sales representative with your specifications.

# MP5




## FEATURES

- Agile & Versatile Mid Wheel Drive Design
- Active Anti-Tippers & Suspension System
- Tight Turning Radius Of 24"
- 350 lbs Load Capacity
- Fully Programmable P & G VSI Controller
- Standard Flip Up Footboard
- Compact Adjustable Height Ergonomically Designed
- Captain Seat
- Rehab Frame Also Available
- Flip Back Arms

SPECIFICATIONS	MP5
OVERALL LENGTH:	42"
OVERALL WIDTH:	25"
WEIGHT CAPACITY:	350 lbs.
SEATING:	18"/20" CAPTAIN SEAT. REHAB FRAME 14" TO 26"
BATTERY SPECS.:	12V 50A x 2 PCS
CONTROLLER TYPE:	P & G FULL PROGRAMMABLE
CHARGER TYPE:	5 AMP, OFF BOARD 120/240 VOLT
MAX SPEED:	8 km/h OR 5 MPH
WEIGHT WITH BATTERY:	215 lbs.
WEIGHT WITHOUT BATTERY:	150 lbs.
DRIVE WHEEL:	12" x 4" SEMI PNEUMATIC
FRONT CASTER:	8" x 2" SEMI PNEUMATIC
BATTERY RANGE:	35 km (22 mi)
MOTOR TYPE:	4600 RPM 450W 32:1
SUSPENSION:	REAR CASTER
TURNING RADIUS:	24"

**MP5**  
POWER WHEELCHAIR



PP-PWCMP5

# HORIZON

NON-TILT COMMODE

STAINLESS STEEL



Power Plus Mobility's **HORIZON** Commode is promised to excel in quality and functionality. With heavy-duty material and ensured high quality assembly, our **HORIZON** is a guarantee that it satisfies your every demand in a commode. Our products will move you in a way that we aim our goals to fulfill any expectations and needs that our customers want in a wheelchair such as our **HORIZON**.

## Horizon 1

### All Features Of Horizon Are The Following :

- Weight Capacity 250 Lbs
- 18" X 18" Easy to Clean Urethane Seat
- Finished Seat Height Including Seat Is 21"
- Seat width By Frame 18"
- Overall Width 20"
- Rear Axle Receiver Ready
- Stainless Steel Frame - Water Resistant
- 20" Back Height
- Black Urethane Easy To Clean Back Rest
- Black Molded Sanitary Cover
- Flip Back Arms
- Black Urethane Easy To Clean Arm Pads
- 6" X 10" X 14" Large Bucket
- 4" X 5" Rust Resistant Nylon Casters (3 Total Locking Casters and 1 Directional Lock Caster)
- Water Proof Headrest

## Upcharge Options

- 14" To 20" Seat Depth
- 21" To 26" Seat Depth
- 14" To 17" Seat Width
- 23" To 26" Seat Width
- Custom Back Height
- Capacity 750 To 1000 Lbs
- Tension Back
- Lock Down Flip Back Arms
- Closed Oval Soft Seat 18" To 22"
- Key Hole Soft Seat 18" To 22"
- 60° Footrest
- 70° Footrest
- 80° Footrest
- 90° Footrest
- Elevating Legrest (Not Water Proof)
- Angle Adjustable Footplates (Not Water Proof)
- One Piece Footboard
- Seat Belt (Plastic)

### Seat Heights

- 24" Wheels X 5" Casters = 21" Seat to Floor Height

### Wheel Self Propelling Wheels:

- 24" Wheels/ 5" Casters

### Warranty: 5 Years On Frame

- 1 Year On Components



PP-MWCHOR

## Horizon 2

### All Features Of Horizon 1 Plus The Following:

- Weight Capacity 500 Lbs
- 18" X 18" Easy to Clean Urethane Seat
- Seat Width By Frame 20"
- Seat Width By Frame 20"
- Overall With 22"
- Heavy Duty Construction

## Horizon 3

### All Features Of Horizon 1 Plus The Following:

- Weight Capacity 750 Lbs
- 21" X 18" Easy to Clean Urethane Seat
- Seat Width By Frame 22"
- Overall With 24"
- Extra Heavy Duty Construction



POWERPLUSMOBILITY.COM

HORIZON  
COMMODE



# HORIZON TILT

TILT COMMODE

STAINLESS STEEL



Power Plus Mobility introduces **HORIZON TILT**. Bringing a commode with far more accessibility and maneuverability than competitors, this product will surely be one of your top choices. Commodes are essentially a product that needs to meet a variety of requirements because of its use. Power Plus Mobility's Commode line will surely fulfill all requirements and needs that are expected from a product of this nature.

## Horizon Tilt 1

### All Features Of Horizon Are The Following :

- Weight Capacity 250 Lbs
- 18" X 18" Easy to Clean Urethane Seat
- Finished Seat Height Including Seat Is 21"
- Seat Width By Frame 18"
- Overall Width 20"
- Rear Axle Receiver Ready
- Silver Powder Coated Finish
- Black Urethane Contour Headrest
- 30° Of Tilt
- Black Urethane Easy to Clean Back Rest
- Black Molded Sanitary Cover
- Flip Back Arms
- Black Urethane Easy To Clean Arm Pads
- 6" X 10" X 14" Large Bucket
- 20" Back Height
- 4" X 5" Rust Resistant Nylon Casters (3 Total Locking Casters And 1 Directional Lock Caster)

## Upcharge Options

- 14" To 20" Seat Depth
- 21" To 26" Seat Depth
- 14" To 17" Seat Width
- 23" To 26" Seat Width
- Custom Back Height
- Capacity 750 To 1000 Lbs
- Tension Back
- Lock Down Flip Back Arms
- Custom Soft Seat 18" To 22" (Closed Oval)
- Custom Soft Seat 18" To 22" (Key Hole)
- 60° Footrest
- 70° Footrest
- 80° Footrest
- 90° Footrest
- Elevating Legrest - Not Water Proof
- Angle Adjustable Footplates - Not Water Proof
- One Piece Footboard
- Seat Belt (Plastic)

### Self Propelling Wheels:

- 24" Wheels
- 22" Wheels
- 20" Wheels
- Rear Wheels Only



### All Features Of Horizon 1 Plus The Following:

- Weight Capacity 500 Lbs
- 18" X 18" Easy to Clean Urethane Seat
- Seat Width By Frame 20"
- Overall Width 22"
- Heavy Duty Construction


## Horizon Tilt 3

### All Features Of Horizon 1 Plus The Following:

- Weight Capacity 750 Lbs
- 21" X 18" Easy to Clean Urethane Seat
- Seat Width By Frame 22"
- Overall With 24"
- Extra Heavy Duty Construction

### Warranty: 5 Years On Frame

- 1 Year On Components



PP-TWCHOR



# CUSTOMIZE ANYTHING

**POWER PLUS MOBILITY** advanced machinery specializes in manufacturing custom parts to customers specification from either a drawing or a sample. The parts we manufacture could be simple or more complex and can be manufactured in small or large quantities.

Contact your local sales representative with your specifications.





Established in December 2002  
located in Ontario, Canada;  
**Power Plus Mobility** strives to  
further enhance the condition  
of life for all mobility  
challenged individuals by  
always making the customer  
our first priority. We  
continuously aspire to  
manufacture and distribute  
improved innovative mobility  
products and seating  
systems. Promoting whatever  
you desire and intended use  
may be, whether it's for an  
active lifestyle or clinical care.  
We promise to provide our  
clients with the best possible  
experience because at Power  
Plus Mobility, we ensure our  
products will truly move you.

## WHEELCHAIR BOOKLET



**POWERPLUSMOBILITY.COM**

208 WILKINSON ROAD | BRAMPTON, ON | CANADA L6T 4M4  
TELEPHONE (905)-614-0333 | TOLL FREE 1-(855)-614-0333  
FAX (905)-614-0435